

TMB-TUNGSTATE METHOD FROM IDA LLEWELLYN-SMITH

1. Incubate tissue in following solution for 10 mins.
10ml 0.1M Phosphate buffer pH 6.0. (see note 1)
0.5ml 1% Ammonium p-tungstate (see note 2)
0.25 ml 0.2% TMB in absolute alcohol (see note 3)
2. Add 10 μ l 3% H₂O₂ and let tissue react for up to 1 hour. (see note 4)
3. Wash in 0.1M phosphate pH 6.0 (2 x 5mins)

This TMB reaction product is stable in everything except osmium so you can mount some sections for LM, dehydrate as normal and coverslip.

TO STABILIZE FOR EM

1. Incubate tissue in following solution for 10 mins
10ml 0.1M Phosphate buffer pH 6.0.
200 μ l 1% Cobalt Chloride
200 μ l 25 mg/ml DAB in water
200 μ l 0.3% H₂O₂
2. Wash in 0.1M phosphate pH 6.0 (3 x 5mins)
3. 1% Osmium Tetroxide in phosphate pH 6.0 (1 hour)
4. Rest of procedure as normal

NOTES

1. Dilute stock 0.2M buffer with H₂O and add HCl to bring to pH 6.0
2. Takes a long time to dissolve (>30 mins), Use magnetic stirrer. Can be stored at **room temperature** for about 1 week.
3. This is 0.01g TMB in 5ml alcohol. (According to protocol this solution can be stored for 2 months at 4°C.)
4. 3% H₂O₂ = 1ml H₂O₂ in 10ml dist H₂O₂. (stock H₂O₂ = 30%)