

Fos immunocytochemistry

Day 1

1. Prepare the **BUFFER**: PBS or TPBS with 1% NGS and 0.3% triton.
2. Add antibody (one 20 µl/tube of 1:20 dilution, in Revco freezer) to 0.5 ml of buffer in a glass tube.
3. Place 15-20 mg of rat liver powder¹ in the glass tube.
4. Incubate in the refrigerator overnight or put glass tube in 37^o C oven for 1 hour followed by 1 hour in 4^oC refrigerator. Agitate gently once or twice during this time.
5. Wash tissue + **control** twice in the buffer with agitation.
6. Remove liquid from wells.
7. Pre-incubate the sections, on the shaker, for at least 60 min in **BLOCKING SOLUTION**. DO NOT WASH AFTER BLOCKING.
8. Rinse glass tube containing the antibody-liver powder solution with 5 ml of buffer and transfer this solution (5.5 ml) in a Centrex centrifuge tube. Centrifuge at 3500 rpm for 10 min. Then add 14.5 ml of buffer (final volume =20 ml).
9. Empty the wells and apply 400 µl of primary antibody per well. Incubate for 48 hours on shaker table at 4^o C.

¹ Cooper Biomedical Rat Liver Powder. Cat#0013-1140.

Day 2

1. Wash with buffer 3 times with agitation for 5-10 min between rinses.
2. Dilute secondary antibody (**bridge**): 1 drop of GAR per 10 ml of buffer.
3. Empty the wells and apply 400 μ l of secondary antibody per well.
Incubate for 60 min at room temperature.
4. Wash 3 times with **NEW** buffer: PBS or TPBS with 0.3% triton (NO SERUM).
5. Dilute standard ABC (Vector PK4000) 30 min. before using it: add 2 drops of solution A and then 2 drops of solution B per 10 ml of **new** buffer. Never mix A and B before adding the buffer. Never mix solutions from different kits.
6. Empty wells and apply 400 μ l of ABC per well, for 60 min with agitation.
7. Wash three times with PBS or TPBS (no triton, no serum).
8. Refrigerate until DAB reaction.