

Thursday, March 19, 1998

## Protocol for Histamine-Elisa Kit from IBL

Before starting check if:

- 1- Is the **assay buffer** prepared? Dilute the assay buffer concentrate 1:5 with dd H<sub>2</sub>O. Store at 2-8°C
- 2- Is the **wash buffer** prepared? Dilute the wash buffer concentrate 1:20 with dd H<sub>2</sub>O.
- 3- Is the **enzyme conjugate** ready? 5µl of enzyme conjugate concentrate in 1.0 ml of assay buffer. Must be made fresh every day.

### I- ACYLATE STANDARDS, CONTROLS, AND SPECIMENS

- 1- Get one 3 ml glass tube per sample and identify
- 2- Fill glass tubes with:  
**50 µl** of U/Z standards (A to F) or **control urine** or **rat urine** in appropriate tube  
**50 µl** of **indicator buffer** in each tube (must stay pink otherwise adjust pH)  
**10 µl** of **acetylation reagent** in each tube
- 3- Incubate a room temp for **30 min.**
- 4- Add 2 ml of **assay buffer** (*diluted*) and vortex mix.

### II- IMMUNOREACTION

- 1- Pipet **50 µl** of acetylated standards, control, and rat urine in appropriate well
- 2- Add **50 µl** of **enzyme conjugate** (*use diluted conjugate!*)
- 3- Add **50 µl** of **antiserum** (green)
- 4- Seal and incubate on shaker for **3 hours**

### III- WASHING

Wash each well 4 times with 250 µl of **wash buffer** (*diluted*)

### IV- TMB reaction

- 1- Prepare TMB substrate solution: 100µl of **TMB substrate concentrate** + 3 ml of **TMB Substrate buffer** and mix.
- 2- Pipet **200µl per well**
- 3- Seal and incubate for **20 min** (urine) on shaker
- 4- Stop by adding 100µl of stop solution into each well.

### V- Read density at 450 nm